

MK12 Polyclonal Antibody

Cat No: HR1AP11671

For research use only

Overview

Product Name	MK12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	MAPK12 ERK6 SAPK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mitogen-activated protein kinase 12(MAPK12) Homo sapiens Activation of members of the mitogen-activated protein kinase family is a major mechanism for transduction of extracellular signals. Stress-activated protein kinases are one subclass of MAP kinases. The protein encoded by this gene functions as a signal transducer during differentiation of myoblasts to myotubes. [provided by RefSeq, Jul 2008],