

MAGG1 Polyclonal Antibody

Cat No: HR1AP11669

For research use only

Overview

Product Name	MAGG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	NDNL2 HCA4 MAGEG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NSE3 homolog, SMC5-SMC6 complex component(NSMCE3) Homo sapiens The protein encoded by this gene is part of the SMC5-6 chromatin reorganizing complex and is a member of the MAGE superfamily. This is an intronless gene. [provided by RefSeq, May 2011],