

MAGD4 Polyclonal Antibody

Cat No: HR1AP11667

For research use only

Overview

Product Name	MAGD4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	MAGED4 KIAA1859 MAGED4A MAGEE1; MAGED4B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAGE family member D4B(MAGED4B) Homo sapiens This gene is a member of the MAGED gene family. It is expressed only in brain and ovary, and some transcript variants of this gene are specifically expressed in glioma cells. This gene is clustered with other MAGED genes on chromosome Xp11. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, May 2011],