

MAGE2 Polyclonal Antibody

Cat No: HR1AP11668

For research use only

Overview

Product Name	MAGE2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	MAGEE2 HCA3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAGE family member E2(MAGEE2) Homo sapiens This gene encodes a member of the E subfamily of MAGE (melanoma antigen-encoding gene) gene family. The gene is intronless and the encoded protein has two of the MAGE domains which are characteristic of MAGE family proteins. [provided by RefSeq, Sep 2011],