

MAGD2 Polyclonal Antibody

Cat No: HR1AP11666

For research use only

Overview

Product Name	MAGD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	MAGED2 BCG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAGE family member D2(MAGED2) Homo sapiens This gene is a member of the MAGED gene family. The MAGED genes are clustered on chromosome Xp11. This gene is located in Xp11.2, a hot spot for X-linked mental retardation (XLMR). This gene may also be involved in several types of cancer, including breast cancer and melanoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013].