

# MAGBI Polyclonal Antibody

Cat No: HR1AP11664

For research use only

## Overview

Product Name	MAGBI Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	MAGEB18
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 MAGE domain.,