

MAGBA Polyclonal Antibody

Cat No: HR1AP11663

For research use only

Overview

Product Name	MAGBA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	MAGEB10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAGE family member B10(MAGEB10) Homo sapiens This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells. [provided by RefSeq, Apr 2010],