

## CUL7 Polyclonal Antibody

Cat No: HR1AP11645

For research use only

### Overview

Product Name	CUL7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	186kDa
GeneID?Human?	CUL7 KIAA0076
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cullin 7(CUL7) Homo sapiens The protein encoded by this gene is a component of an E3 ubiquitin-protein ligase complex. The encoded protein interacts with TP53, CUL9, and FBXW8 proteins. Defects in this gene are a cause of 3M syndrome type 1 (3M1). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009],</p>