

CYLC1 Polyclonal Antibody

Cat No: HR1AP11646

For research use only

Overview

Product Name	CYLC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	71kDa
GeneID?Human?	CYLC1 CYL CYL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cylicin 1(CYLC1) Homo sapiens This gene encodes a sperm head cytoskeletal protein. The encoded protein is associated with the calyx of spermatozoa and spermatids. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012],