



## CUL4A Polyclonal Antibody

Cat No: HR1AP11643

For research use only

### Overview

Product Name	CUL4A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	CUL4A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cullin 4A(CUL4A) Homo sapiens CUL4A is the ubiquitin ligase component of a multimeric complex involved in the degradation of DNA damage-response proteins (Liu et al., 2009 [PubMed 19481525]).[supplied by OMIM, Oct 2009],