

CENPO Polyclonal Antibody

Cat No: HR1AP11638

For research use only

Overview

Product Name	CENPO Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	CENPO ICEN36 MCM21R
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	centromere protein O(CENPO) Homo sapiens This gene encodes a component of the interphase centromere complex. The encoded protein is localized to the centromere throughout the cell cycle and is required for bipolar spindle assembly, chromosome segregation and checkpoint signaling during mitosis. Alternatively spliced transcript variants encoding multiple protein isoforms have been observed for this gene. [provided by RefSeq, Dec 2010],