

CDK12 Polyclonal Antibody

Cat No: HR1AP11629

For research use only

Overview

Product Name	CDK12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	163kDa
GeneID?Human?	CDK12 CRK7 CRKRS KIAA0904
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the protein kinase superfamily.,