

CDK13 Polyclonal Antibody

Cat No: HR1AP11627

For research use only

Overview

Product Name	CDK13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	166kDa
GeneID?Human?	CDK13 CDC2L CDC2L5 CHED KIAA1791
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cyclin dependent kinase 13(CDK13) Homo sapiens The protein encoded by this gene is a member of the cyclin-dependent serine/threonine protein kinase family. Members of this family are well known for their essential roles as master switches in cell cycle control. The exact function of this protein has not yet been determined, but it may play a role in mRNA processing and may be involved in regulation of hematopoiesis. Alternatively spliced transcript variants have been described.[provided by RefSeq, Dec 2009],</p>