



CCNY Polyclonal Antibody

Cat No: HR1AP11624

For research use only

Overview

Product Name	CCNY Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	CCNY C10orf9 CBCP1 CFP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cyclin Y(CCNY) Homo sapiens Cyclins, such as CCNY, control cell division cycles and regulate cyclin-dependent kinases (e.g., CDC2; MIM 116940) (Li et al., 2009 [PubMed 18060517]).[supplied by OMIM, May 2009],