

CC50A Polyclonal Antibody

Cat No: HR1AP11615

For research use only

Overview

Product Name	CC50A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	TMEM30A C6orf67 CDC50A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the CDC50/LEM3 family.,