

ASPP2 Polyclonal Antibody

Cat No: HR1AP11611

For research use only

Overview

Product Name	ASPP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	124kDa
GeneID?Human?	TP53BP2 ASPP2 BBP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tumor protein p53 binding protein 2(TP53BP2) Homo sapiens This gene encodes a member of the ASPP (apoptosis-stimulating protein of p53) family of p53 interacting proteins. The protein contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. It is localized to the perinuclear region of the cytoplasm, and regulates apoptosis and cell growth through interactions with other regulatory molecules including members of the p53 family. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],