

ARAID Polyclonal Antibody

Cat No: HR1AP11610

For research use only

Overview

Product Name	ARAID Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	ATRAID APR3 C2orf28 HSPC013 UNQ214/PRO240
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	all-trans retinoic acid induced differentiation factor(ATRAID) Homo sapiens This gene is thought to be involved in apoptosis, and may also be involved in hematopoietic development and differentiation. The use of alternative splice sites and promoters result in multiple transcript variants encoding different isoforms.[provided by RefSeq, Dec 2009],