

ADCY2 Polyclonal Antibody

Cat No: HR1AP11609

For research use only

Overview

Product Name	ADCY2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	120kDa
GeneID?Human?	ADCY2 KIAA1060
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	adenylate cyclase 2(ADCY2) Homo sapiens This gene encodes a member of the family of adenylate cyclases, which are membrane-associated enzymes that catalyze the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This enzyme is insensitive to Ca(2+)/calmodulin, and is stimulated by the G protein beta and gamma subunit complex. [provided by RefSeq, Jul 2008],