

ZA2G Polyclonal Antibody

Cat No: HR1AP11608

For research use only

Overview

Product Name	ZA2G Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	AZGP1 ZAG ZNGP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acids.,similarity:Belongs to the MHC class I family.,similarity:Contains 1 Ig-like C1-type (immunoglobulin-like) domain.,subunit:Interacts with PIP.,tissue specificity:Blood plasma, seminal plasma, urine, saliva, sweat, epithelial cells of various human glands, liver.,