

S39A6 Polyclonal Antibody

Cat No: HR1AP11605

For research use only

Overview

Product Name	S39A6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	SLC39A6 LIV1 ZIP6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 39 member 6(SLC39A6) Homo sapiens Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008],</p>