

ZNT3 Polyclonal Antibody

Cat No: HR1AP11604

For research use only

Overview

Product Name	ZNT3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	SLC30A3 ZNT3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Involved in accumulation of zinc in synaptic vesicles.,similarity:Belongs to the cation diffusion facilitator (CDF) transporter (TC 2.A.4) family. SLC30A subfamily.,