

## ZHX3 Polyclonal Antibody

Cat No: HR1AP11603

For research use only

### Overview

Product Name	ZHX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	ZHX3 KIAA0395 TIX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>zinc fingers and homeoboxes 3(ZHX3) Homo sapiens This gene encodes a member of the zinc fingers and homeoboxes (ZHX) gene family. The encoded protein contains two C2H2-type zinc fingers and five homeodomains and forms a dimer with itself or with zinc fingers and homeoboxes family member 1. In the nucleus, the dimerized protein interacts with the A subunit of the ubiquitous transcription factor nuclear factor-Y and may function as a transcriptional repressor. [provided by RefSeq, Jul 2008].</p>