

ZN462 Polyclonal Antibody

Cat No: HR1AP11601

For research use only

Overview

Product Name	ZN462 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	275kDa
GeneID?Human?	ZNF462 KIAA1803
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 27 C2H2-type zinc fingers.,