

## ZN335 Polyclonal Antibody

Cat No: HR1AP11600

For research use only

### Overview

Product Name	ZN335 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	147kDa
GeneID?Human?	ZNF335
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger protein 335(ZNF335) Homo sapiens The protein encoded by this gene enhances transcriptional activation by ligand-bound nuclear hormone receptors. However, it does this not by direct interaction with the receptor, but by direct interaction with the nuclear hormone receptor transcriptional coactivator NRC. The encoded protein may function by altering local chromatin structure. [provided by RefSeq, Jul 2008],