

ZN238 Polyclonal Antibody

Cat No: HR1AP11598

For research use only

Overview

Product Name	ZN238 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	ZNF238 RP58 TAZ1 ZBTB18
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger and BTB domain containing 18(ZBTB18) Homo sapiens This gene encodes a C2H2-type zinc finger protein which acts a transcriptional repressor of genes involved in neuronal development. The encoded protein recognizes a specific sequence motif and recruits components of chromatin to target genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013],