

ZFH2 Polyclonal Antibody

Cat No: HR1AP11595

For research use only

Overview

Product Name	ZFH2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	282kDa
GeneID?Human?	ZFH2 KIAA1056 KIAA1762 ZNF409
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in transcriptional regulation.,similarity:Contains 3 homeobox DNA-binding domains.,similarity:Contains 6 C2H2-type zinc fingers.,