

TUT4 Polyclonal Antibody

Cat No: HR1AP11593

For research use only

Overview

Product Name	TUT4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	180kDa
GeneID?Human?	ZCCHC11 KIAA0191 TUT4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger CCHC-type containing 11(ZCCHC11) Homo sapiens ZCCHC11 is an RNA uridylyltransferase (EC 2.7.7.52) that uses UTP to add uridines to the 3-prime end of substrate RNA molecules (Jones et al., 2009 [PubMed 19701194]).[supplied by OMIM, Jan 2011],