

ZBTB2 Polyclonal Antibody

Cat No: HR1AP11591

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | ZBTB2 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 56kDa |
| GeneID?Human? | ZBTB2 KIAA1483 ZNF437 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers., |