

WDR44 Polyclonal Antibody

Cat No: HR1AP11588

For research use only

Overview

Product Name	WDR44 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	WDR44
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	WD repeat domain 44(WDR44) Homo sapiens This gene encodes a protein that interacts with the small GTPase rab11. A similar protein in rat binds the GTP-containing active form of rab11. This protein may play a role in endosome recycling. Alternate splicing results in multiple transcript variants.[provided by RefSeq, May 2010],