

CAC1S Polyclonal Antibody

Cat No: HR1AP11583

For research use only

Overview

Product Name	CAC1S Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	206kDa
GeneID?Human?	CACNA1S CACH1 CACN1 CACNL1A3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calcium voltage-gated channel subunit alpha1 S(CACNA1S) Homo sapiens This gene encodes one of the five subunits of the slowly inactivating L-type voltage-dependent calcium channel in skeletal muscle cells. Mutations in this gene have been associated with hypokalemic periodic paralysis, thyrotoxic periodic paralysis and malignant hyperthermia susceptibility. [provided by RefSeq, Jul 2008],