

## SFT2A Polyclonal Antibody

Cat No: HR1AP11579

For research use only

### Overview

Product Name	SFT2A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	SFT2D1 C6orf83
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,