

VAX2 Polyclonal Antibody

Cat No: HR1AP11576

For research use only

Overview

Product Name	VAX2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	VAX2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ventral anterior homeobox 2(VAX2) Homo sapiens This gene encodes a homeobox protein and is almost exclusively expressed in the ventral portion of the retina during development. In mouse studies, this gene was found to be required for the correct formation of the optic fissure and other aspects of retinal development. [provided by RefSeq, Sep 2008],