

VANG2 Polyclonal Antibody

Cat No: HR1AP11574

For research use only

Overview

Product Name	VANG2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	VANGL2 KIAA1215 STB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	VANGL planar cell polarity protein 2(VANGL2) Homo sapiens The protein encoded by this gene is a membrane protein involved in the regulation of planar cell polarity, especially in the stereociliary bundles of the cochlea. The encoded protein transmits directional signals to individual cells or groups of cells in epithelial sheets. This protein is also involved in the development of the neural plate. [provided by RefSeq, Sep 2011],