

SOGA2 Polyclonal Antibody

Cat No: HR1AP11562

For research use only

Overview

Product Name	SOGA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	209kDa
GeneID?Human?	SOGA2 CCDC165 KIAA0802
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Isoform 4 is phosphorylated upon DNA damage, probably by ATM or ATR.,