

RPRD2 Polyclonal Antibody

Cat No: HR1AP11561

For research use only

Overview

Product Name	RPRD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	160kDa
GeneID?Human?	RPRD2 KIAA0460 HSPC099
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 CID domain.,