

UBXN4 Polyclonal Antibody

Cat No: HR1AP11559

For research use only

Overview

Product Name	UBXN4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	UBXN4 KIAA0242 UBXD2 UBXDC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	UBX domain protein 4(UBXN4) Homo sapiens UBXD2 is an integral membrane protein of the endoplasmic reticulum (ER) that binds valosin-containing protein (VCP; MIM 601023) and promotes ER-associated protein degradation (ERAD) (Liang et al., 2006 [PubMed 16968747]).[supplied by OMIM, Mar 2008],