

UBE2Z Polyclonal Antibody

Cat No: HR1AP11557

For research use only

Overview

Product Name	UBE2Z Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	UBE2Z HOYS7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ubiquitin conjugating enzyme E2 Z(UBE2Z) Homo sapiens This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis. [provided by RefSeq, Feb 2012],