

UB2R2 Polyclonal Antibody

Cat No: HR1AP11554

For research use only

Overview

Product Name	UB2R2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	UBE2R2 CDC34B UBC3B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ubiquitin conjugating enzyme E2 R2(UBE2R2) Homo sapiens Protein kinase CK2 is a ubiquitous and pleiotropic Ser/Thr protein kinase involved in cell growth and transformation. This gene encodes a protein similar to the E2 ubiquitin conjugating enzyme UBC3/CDC34. Studies suggest that CK2-dependent phosphorylation of this ubiquitin-conjugating enzyme functions by regulating beta-TrCP substrate recognition and induces its interaction with beta-TrCP, enhancing beta-catenin degradation. [provided by RefSeq, Jul 2008],</p>