

UB2Q1 Polyclonal Antibody

Cat No: HR1AP11552

For research use only

Overview

Product Name	UB2Q1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	UBE2Q1 NICE5 UBE2Q PRO3094
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ubiquitin conjugating enzyme E2 Q1(UBE2Q1) Homo sapiens The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s), and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein is 98% identical to the mouse counterpart. [provided by RefSeq, Jul 2008].</p>