

## UBAP2 Polyclonal Antibody

Cat No: HR1AP11542

For research use only

### Overview

Product Name	UBAP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	123kDa
GeneID?Human?	UBAP2 KIAA1491
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ubiquitin associated protein 2(UBAP2) Homo sapiens The protein encoded by this gene contains a UBA (ubiquitin associated) domain, which is characteristic of proteins that function in the ubiquitination pathway. This gene may show increased expression in the adrenal gland and lymphatic tissues. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013],