

# UB2L6 Polyclonal Antibody

Cat No: HR1AP11539

For research use only

## Overview

Product Name	UB2L6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	UBE2L6 UBCH8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ubiquitin conjugating enzyme E2 L6(UBE2L6) Homo sapiens The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s) and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is highly similar in primary structure to the enzyme encoded by the UBE2L3 gene. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2011],</p>