

UBP6 Polyclonal Antibody

Cat No: HR1AP11535

For research use only

Overview

Product Name	UBP6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	154kDa
GeneID?Human?	USP6 HRP1 TRE2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,disease:Involved in tumorigenesis.,function:Has an ATP-independent isopeptidase activity, cleaving at the C-terminus of the ubiquitin moiety. In vitro, isoform 2, but not isoform 3, shows deubiquitinating activity.,miscellaneous:The USP6 gene only exists in the primate lineage.,similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 Rab-GAP TBC domain.,tissue specificity:Testis specific. Expressed in various cancer cell lines.,</p>