

# UBP37 Polyclonal Antibody

Cat No: HR1AP11531

For research use only

## Overview

Product Name	UBP37 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	107kDa
GeneID?Human?	USP37 KIAA1594
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,similarity:Contains 3 UIM (ubiquitin-interacting motif) repeats.,tissue specificity:Expressed in brain and prostate.,