

UBP31 Polyclonal Antibody

Cat No: HR1AP11528

For research use only

Overview

Product Name	UBP31 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	148kDa
GeneID?Human?	USP31 KIAA1203
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,function:May recognize and hydrolyze the peptide bond at the C-terminal Gly of ubiquitin. Involved in the processing of poly-ubiquitin precursors as well as that of ubiquitinated proteins.,similarity:Belongs to the peptidase C19 family.,tissue specificity:Widely expressed.,</p>