



## UBP3 Polyclonal Antibody

Cat No: HR1AP11527

For research use only

### Overview

Product Name	UBP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	USP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,similarity:Contains 1 UBP-type zinc finger.,tissue specificity:Expressed in all tissues examined, with strongest expression in pancreas.,