

## UBP26 Polyclonal Antibody

Cat No: HR1AP11525

For research use only

### Overview

Product Name	UBP26 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	USP26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ubiquitin specific peptidase 26(USP26) Homo sapiens This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility. [provided by RefSeq, Jul 2008],