

UBP1 Polyclonal Antibody

Cat No: HR1AP11520

For research use only

Overview

Product Name	UBP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	USP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ubiquitin specific peptidase 1(USP1) Homo sapiens This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases that is a deubiquitinating enzyme (DUB) with His and Cys domains. This protein is located in the cytoplasm and cleaves the ubiquitin moiety from ubiquitin-fused precursors and ubiquitinated proteins. The protein specifically deubiquitinates a protein in the Fanconi anemia (FA) DNA repair pathway. Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008].</p>