

UBQL2 Polyclonal Antibody

Cat No: HR1AP11519

For research use only

Overview

Product Name	UBQL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	UBQLN2 N4BP4 PLIC2 HRIHFB2157
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ubiquilin 2(UBQLN2) Homo sapiens This gene encodes an ubiquitin-like protein (ubiquilin) that shares high degree of similarity with related products in yeast, rat and frog. Ubiquilins contain a N-terminal ubiquitin-like domain and a C-terminal ubiquitin-associated domain. They physically associate with both proteasomes and ubiquitin ligases; and thus, are thought to functionally link the ubiquitination machinery to the proteasome to affect in vivo protein degradation. This ubiquilin has also been shown to bind the ATPase domain of the Hsp70-like Stch protein. [provided by RefSeq, Oct 2009],</p>