

LSM2 Polyclonal Antibody

Cat No: HR1AP11518

For research use only

Overview

Product Name	LSM2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	LSM2 C6orf28 G7B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>LSM2 homolog, U6 small nuclear RNA and mRNA degradation associated(LSM2) Homo sapiens This gene encodes a member of the LSm family of RNA-binding proteins. LSm proteins form stable heteromers that bind specifically to the 3'-terminal oligo(U) tract of U6 snRNA and may play a role in pre-mRNA splicing by mediating U4/U6 snRNP formation. Pseudogenes of this gene are located on the short arm of chromosomes 6 and 19. [provided by RefSeq, Nov 2011],</p>