

TNNC2 Polyclonal Antibody

Cat No: HR1AP11509

For research use only

Overview

Product Name	TNNC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	TNNC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	troponin C2, fast skeletal type(TNNC2) Homo sapiens Troponin (Tn), a key protein complex in the regulation of striated muscle contraction, is composed of 3 subunits. The Tn-I subunit inhibits actomyosin ATPase, the Tn-T subunit binds tropomyosin and Tn-C, while the Tn-C subunit binds calcium and overcomes the inhibitory action of the troponin complex on actin filaments. The protein encoded by this gene is the Tn-C subunit. [provided by RefSeq, Jul 2008],